

Ophthalmic

INSTRUMENT
CATALOG



Your Care, Our Focus

Contents

SMILE Products	1
Featured Products	2
Speculums	6
Cannulas	10
Forceps	14
Scissors	32
Needle Holder	39
Caliper and Marker	42
Rongeur/Punch/Trephine	47
Hooks/Lacrimonal/Retractor	51
Phaco Chopper/Manipulator/Spatula	55
I/A Handpiece and Flute Needle	64
Sterilizing Cases	67
Trypan Blue	70
Knives	72
Surgical Package	75

F12-001

Fixation Ring

滾滾長江東逝水



- Diameter 13mm
- Overall length: 107mm, 4.2 inches
- Stainless Steel



F12-002

Lenticule Remove Forceps


- Serrated Tips, 23G
- Curved
- Overall length: 129mm, 5 inches
- Titanium



F12-003

Smile Dissector

- Smile Double-Ended Dissector
- Overall length: 127mm, 5 inches
- Stainless Steel



F12-006

Smile Lenticule Double-Ended Splitter



- One end aims at helping surgeons to extract & remove the lenticule
- Overall length: 127mm, 5 inches
- Stainless Steel



F3-302T

Inamura Capsulorhexis Forceps



- 10mm curved shafts, clear marking on tips
- Cross action round handle
- Design down to 1.5mm incision, Small blunt tips insure easy entry
- Overall length: 115mm, 4.5 inches
- Titanium



F3-300T

Inamura Capsulorhexis Forceps


- 11mm curved shafts, clear marking on tips
- Cross action flat handle
- Design down to 1.8mm incision, Small blunt tips insure easy entry
- Overall length: 115mm, 4.5 inches
- Titanium



F9-605

Polisher Double-ended 27g

- 27g
- Designed for removing LECs from the inner surface of the capsule
- The circular ends are ideal for cleaning the equator of the capsular bag prior to lens insertion
- Provides a left angled and right angled polisher on one handle for surgeon convenience
- 1mm diameter both sides
- Overall length: 134mm, 5.3 inches
- Stainless Steel



F3-288T

Micro Incision Forceps


- 12mm curved shafts
- Design down to 1.8mm incision, Small blunt tips insure easy entry
- Overall length: 113mm, 4.4 inches
- Titanium



F3-227T

Micro Incision Forceps


- 12mm straight shafts
- Design down to 1.8mm incision, Small blunt tips insure easy entry
- Overall length: 113mm, 4.4 inches
- Titanium



F3-550

Femto Prechopper


- Overall length: 120mm, 4.7 inches
- Stainless Steel



F3-549

Akahoshi Prechopper


- 23mm straight shafts with crescent tips
- Overall length: 120mm, 4.7 inches
- Stainless Steel
- F3-549T (Titanium)



F3-599

Akahoshi Prechopper


- 22mm angled shafts with sharp point tips
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F3-688T

Artisan Lens Implantation Forceps

- 7mm angled shafts
- Overall length: 124mm, 4.9 inches
- Titanium



F3-788T

Trocar Ports Grasping Forceps


- 23 gauge
- For Easy Removal of Trocar Port From Sclera
- Overall length: 113mm, 4.4 inches
- Titanium



F3-802T-R

Artiflex Lens Implantation Forceps Right



- Overall length: 107mm, 4.2 inches
- Titanium



F3-802T-L

Artiflex Lens Implantation Forceps Left

- Overall length: 107mm, 4.2 inches
- Titanium



F4-361

Dalk Scissors

- 7mm cut length, from mid-screw to tip 10.5mm
- Used to remove the superficial layers of corneal stroma
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F4-361T (Titanium)



F6-036

Toric Axis Marker


- Overall length: 100mm, 3.9 inches
- Stainless Steel



F6-037

Toric Reference Ring Marker


- Overall length: 105mm, 4.1 inches
- Stainless Steel



F6-038

Corneal Marker



- 12 line
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F9-604

Femto / Smile Double-ended Lamellar

- Dissector allows surgeries to open femto
- incisions prepared by the laser
- Overall length: 130mm, 5.1 inches
- Stainless Steel



F3-163

Donnor Membrane Holder / Dmek

- Adapt for DMEK surgery with broad grasping
- the endothelium membrane of the donor
- Overall length: 100mm, 3.9 inches
- Stainless Steel



F9-759T

Lamella Spear

- 7mm diameter
- Overall length: 122mm, 4.8 inches
- Titanium




F9-907T

Lamella Spear

- 9mm diameter
- Overall length: 122mm, 4.8 inches
- Titanium




Wire
Adjustable
Locking
with Aspiration





F1-001 Wire Speculum

- Closed 15mm blades
- Overall length: 50mm, 2.0 inches
- Single piece construction
- Stainless Steel
- F1-007 (Pediatric, 11mm blades)
- F1-009 (New born, 7mm blades)





F1-003 Wire Speculum

- Closed 15mm blades
- Overall length: 50mm, 2.0 inches
- Single piece construction
- Stainless Steel




F1-005T Wire Speculum

- Open 15mm blades
- Overall length: 50mm, 2.0 inches
- Single piece construction
- Titanium
- F1-005 (Stainless Steel)




F1-011T Barraquer Eye Speculum V-shape

- Open 15mm blades
- V-shape
- Overall length: 75mm, 3.0 inches
- Adjustable with thumb screw
- Titanium
- F1-011 (Stainless Steel)




F1-013 Speculum(New born)

- Closed 3.5mm with 6.5mm blades
- Overall length: 45mm, 1.8 inches
- Stainless Steel
- F1-013-L (Left)
- F1-013-R (Right)




F1-021 Speculum (Translational)

- Closed 15mm blades
- 0-36mm adjustable
- Overall length: 45mm, 1.8 inches
- Stainless Steel




F1-023 Cook Speculum (Translational)

- Solid 15mm blades
- Overall length: 45mm, 1.8 inches
- Adjustable with thumb screw
- Stainless Steel
- F1-027 (11mm blades)
- F1-031 (8mm blades)




F1-026 Lieberman Speculum V-shape

- Open 15mm blades
- V-shape
- Overall length: 85mm, 3.3 inches
- Adjustable with thumb screw
- Stainless Steel




F1-032 Barrette Speculum

- Open 14mm blades
- Overall length: 60mm, 2.4 inches
- Adjustable with thumb screw
- Stainless Steel




F1-035T Kershner Eye Speculum

- Solid 15mm blades
- Overall length: 65mm, 2.6 inches
- Adjustable with thumb screw
- Titanium
- F1-015T (13mm blades)
- F1-025T (Pediatric, 8mm blades)




F1-041 Wire Eye Speculum

- Solid 15mm blades
- Overall length: 50mm, 2.0 inches
- Single piece construction
- Stainless Steel
- F1-043 (10mm blades)
- F1-045 (4mm blades)





F1-058 Eye Speculum

- Open 15mm blades
- Overall length: 75mm, 3.0 inches
- Adjustable with thumb screw
- Suitable for femtosecond laser
- Stainless Steel





F1-061 Barraquer Eye Speculum

- Closed 15mm blades
- Overall length: 75mm, 3.0 inches
- Adjustable with thumb screw
- Stainless Steel
- F1-061T (Titanium)




F1-063 Barraquer Eye Speculum

- Closed 12mm blades
- Blade spread 20mm
- Overall length: 65mm, 2.6 inches
- Adjustable with thumb screw
- Stainless Steel
- F1-063T (Titanium)




F1-065 Barraquer Eye Speculum K-shape

- Open 15mm blades
- Blade spread 20mm
- K-shape
- Overall length: 75mm, 3.0 inches
- Adjustable with thumb screw
- Stainless Steel
- F1-065T (Titanium)




F1-067 Lieberman Eye Speculum Aspirating

- Open 15mm blades
- 4 holes inside blade
- Overall length: 75mm, 3.0 inches, 100mm tube
- Adjustable with thumb screw
- Stainless Steel




F1-068 Bieberman Eye Speculum V-shape

- Open 13mm blades
- V-shape
- 6 holes inside blade
- Overall length: 85mm, 3.3inches,
- 100mm tube
- Adjustable with thumb screw
- Stainless Steel




F1-091 Eye Speculum (Pediatric)

- Solid 11mm blades
- Overall length: 40mm, 1.6 inches
- Adjustable with thumb screw
- Stainless Steel
- F1-093 (4mm blades)




F1-102T Ring Holding Adjustable Speculum

- Closed 15mm blades
- Overall length: 75mm, 3.0 inches
- Adjustable with thumb screw
- Titanium



Cannulas

Air Injection
Anesthesia
Anterior Chamber Maintainers
Infusion
Irrigating Aspirating
Lacrimal



F2-001 Irrigation Needle Straight 25g

- Straight shafts
- 25 gauge tube, 5#
- Overall length: 40mm, 1.6 inches
- Stainless Steel



F2-003 Irrigation Needle Curved 25g

- Curved with angled 90 degree shafts
- 25 gauge tube
- Overall length: 33mm, 1.3 inches
- Stainless Steel



F2-005 Irrigation Needle Angled 25g

- Angled with 135 degree shafts
- 25 gauge tube
- Overall length: 33mm, 1.3 inches
- Stainless Steel



F2-011 Irrigation Needle Angled 22g

- Angled 45 degree shafts with 10mm from tip to bend
- 22 gauge tube
- 5 holes in the tip
- Overall length: 42mm, 1.7 inches
- Stainless Steel



F2-031 Anterior Chamber Maintainer 1 Hole 18g

- 0.9mm shaft with 10mm from tip to bend
- 18 gauge tube
- Uses at least 30cm of silastic tubing
- Overall length: 10mm, 0.4 inches
- Stainless Steel



F2-033 Anterior Chamber Maintainer 3 Holes 18g

- 0.9mm shaft with 10mm from tip to bend
- 18 gauge tube
- 3 holes in the tip
- Overall length: 10mm, 0.4 inches
- Stainless Steel



F2-041 Cystotome Needle 30g

- Angled 60 degree down
- 30 gauge tube
- Overall length: 12mm, 0.5 inches
- Stainless Steel



F2-043 Anesthesia Needle 25g

- Beveled tip
- 25 gauge tube
- Retrobular atkinson
- Overall length: 67mm, 2.6 inches
- Stainless Steel



F2-051 Infusion Cannula for Water Infusion 16g

- 2.5mm beveled tip and luer lock adaptor
- 16 gauge
- Use 160mm silica gel tube
- Overall length: 8.1mm, 0.3 inches
- Stainless Steel
- F2-053 (3mm)
- F2-055 (4mm)
- F2-057 (5mm)
- F2-059 (6mm)



F2-061 Infusion Cannula for Silicon Oil Infusion

- 3mm beveled tip and luer lock adaptor
- Overall length: 8.1mm, 0.3 inches
- Stainless Steel
- F2-063 (4mm)
- F2-065 (5mm)
- F2-067 (6mm)



F2-073 Jensen Capsule Polisher 23g

- 1mm smooth bulbous tip with sandblasted
- 23 gauge tube
- 0.6mm diameter
- Overall length: 40mm, 1.6 inches
- Stainless Steel



F2-081 Infusion Cannula Angled 18g

- Angled shaft 12mm from tip to bend
- Curved blunt tip
- 18 gauge tube
- Overall length: 39mm, 1.5 inches
- Stainless Steel




F2-083 Infusion Cannula Angled 20g

- Angled 150 degree shaft, 14mm from tip to bend
- Flettend blunt tip
- 20 gauge tube
- Overall length: 39mm, 1.5 inches
- Stainless Steel



F2-091 Infusion Lens Loop 18g

- Thin wall shaft with 4mm wide by 6mm long loop
- 18 gauge tube
- Overall length: 40mm, 1.6 inches
- Stainless Steel
- F2-093 (5mm x 7mm)
- F2-095 (3mm x 8mm)
- F2-097 (3mm x 12mm)



F2-111 Irrigation / Aspiration Cannula 23g

- 0.25mm aspiration
- 23 gauge tube
- Overall length: 32mm, 1.3 inches
- Stainless Steel
- F2-113 (0.30mm, 22 gauge)
- F2-115 (0.35mm, 21 gauge)
- F2-117 (0.40mm, 20 gauge)



F2-121 Simcoe I/A Cannula 22g

- Aspiration angled 90 degree left 2mm flat tip
- 22 gauge tube, 7#
- Overall length: 35mm, 1.4 inches
- Stainless Steel
- F2-121-L (Left)
- F2-121-R (Right)




F2-125 Simcoe I/A Cannula 21g

- Aspiration angled 90 degree left 2mm flat tip
- 21 gauge tube
- Overall length: 35mm, 1.4 inches
- Stainless Steel
- F2-125-L (Left)
- F2-125-R (Right)





F2-129 Lacrimal Cannulas 20g

- 20 gauge
- 0.60mm straight shaft inside
- Bowman lacrimal probe mallerble
- Overall length: 84mm, 3.3 inches
- Stainless Steel



F2-131 Lacrimal Cannulas 25g / 23g / 22g / 21g / 20g

- 25 gauge
- 0.20mm aspiration port
- Overall length: 72mm, 2.8 inches
- Stainless Steel
- F2-133 (0.25mm, 23 gauge)
- F2-135 (0.30mm, 22 gauge)
- F2-137 (0.35mm, 21 gauge)
- F2-139 (0.40mm, 20 gauge)



F2-141 Double-ended Lacrimal Probes

- Double-ended
- 25 gauge and 23 gauge
- Overall length: 134mm, 5.3 inches
- Stainless Steel
- F2-143 (22-21 gauge)
- F2-145 (20-19 gauge)
- F2-147 (20-19 gauge, with hole)
- F2-149 (18-17 gauge)



F2-141T Double-ended Lacrimal Probes (5-6)

- Double-ended
- 25 gauge and 23 gauge
- Overall length: 134mm, 5.3 inches
- Titanium
- F2-143T (7-8 gauge)
- F2-145T (9-10 gauge)
- F2-147T (9-10 gauge, with hole)
- F2-149T (11-12 gauge)



F2-161 Irrigation / Aspiration Cannula 22g

- Aspiration angled right 2mm flat tip
- 22 gauge tube
- Overall length: 35mm, 1.4 inches
- Stainless Steel
- F2-161-L (Left)
- F2-161-R (Right)



F2-165 Irrigation / Aspiration Cannula 21g


- 21 gauge
- Aspiration angled right 2mm flat tip
- Stainless Steel
- Overall length: 35mm, 1.4 inches
- Stainless Steel
- F2-165-L (Left)
- F2-165-R (Right)



Forceps


Suturing
Capsulorhexis
Tying
Retinal
Tissue
Utility
Iris
Fixation
Prechopper
Lens Implantation
Lasik Lasek Flap
Strabismus Muscle





90° Tip

- 0.12mm
- 0.3mm
- 0.5mm



45° Tip

- 0.12mm
- 0.3mm
- 0.5mm

F3-003-90 Bonn Suturing Forceps 0.12mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-003T-90 (Titanium)

F3-003-45 Bonn Suturing Forceps 0.12mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-003T-45 (Titanium)

F3-014-90 Bonn Suturing Forceps 0.3mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-014T-90 (Titanium)

F3-014-45 Bonn Suturing Forceps 0.3mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-014T-45 (Titanium)

F3-016-90 Bonn Suturing Forceps 0.5mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-016T-90 (Titanium)

F3-016-45 Bonn Suturing Forceps 0.5mm Straight 106mm

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-016T-45 (Titanium)

F3-019-45 Bonn Suturing Forceps 0.12mm Straight 85mm


- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 85mm, 3.3 inches
- Stainless Steel
- F3-019T-45 (Titanium)

F3-021-45 Bonn Suturing Forceps 0.3mm Straight 85mm

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 85mm, 3.3 inches
- Stainless Steel
- F3-021T-45 (Titanium)


F3-023-45 Bonn Suturing Forceps 0.5mm Straight 85mm

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 85mm, 3.3 inches
- Stainless Steel
- F3-023T-45 (Titanium)



90° Tip

- 0.12mm
- 0.3mm
- 0.5mm



45° Tip

- 0.12mm
- 0.3mm
- 0.5mm

F3-007-90 **Castroviejo Suturing Forceps 0.12mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-007T-90 (Titanium)

F3-007-45 **Castroviejo Suturing Forceps 0.12mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-007T-45 (Titanium)

F3-256-90 **Castroviejo Suturing Forceps 0.3mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-256T-90 (Titanium)

F3-256-45 **Castroviejo Suturing Forceps 0.3mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-256T-45 (Titanium)

F3-255-90 **Castroviejo Suturing Forceps 0.5mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-255T-90 (Titanium)

F3-255-45 **Castroviejo Suturing Forceps 0.5mm Straight 110mm**

- Straight shafts, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-255T-45 (Titanium)

F3-028-45 **Castroviejo Suturing Forceps 0.12mm Straight 85mm**

- Straight shafts, 1 x 2 Teeth
- 0.12mm, with 5mm platform
- Overall length: 85mm, 3.3 inches
- Stainless Steel
- F3-028T-45 (Titanium)

F3-030-45 **Castroviejo Suturing Forceps 0.3mm Straight 85mm**

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 5mm platform
- Overall length: 85mm, 3.3 inches
- Stainless Steel
- F3-030T-45 (Titanium)

F3-032-45 **Castroviejo Suturing Forceps 0.5mm Straight 85mm**

- Straight shafts, 1 x 2 Teeth
- 0.3mm, with 5mm platform
- Overall length 85mm, 3.3 inches
- Stainless Steel
- F3-032T-45 (Titanium)



90° Tip

- 0.12mm
- 0.3mm
- 0.5mm



45° Tip

- 0.12mm
- 0.3mm
- 0.5mm

F3-201-90 Bonn Suturing Forceps 0.12mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-201T-90 (Titanium)

F3-201-45 Bonn Suturing Forceps 0.12mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-201T-45 (Titanium)

F3-208-90 Bonn Suturing Forceps 0.3mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-208T-90 (Titanium)

F3-208-45 Bonn Suturing Forceps 0.3mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-208T-45 (Titanium)

F3-210-90 Bonn Suturing Forceps 0.5mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-210T-90 (Titanium)

F3-210-45 Bonn Suturing Forceps 0.5mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-210T-45 (Titanium)



90° Tip

- 0.12mm
- 0.3mm
- 0.5mm



45° Tip

- 0.12mm
- 0.3mm
- 0.5mm

F3-009-90 Castroviejo Suturing Forceps 0.12mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-009T-90 (Titanium)

F3-009-45 Castroviejo Suturing Forceps 0.12mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.12mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-009T-45 (Titanium)

F3-006-90 Castroviejo Suturing Forceps 0.3mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-006T-90 (Titanium)

F3-006-45 Castroviejo Suturing Forceps 0.3mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.3mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-006T-45 (Titanium)

F3-008-90 Castroviejo Suturing Forceps 0.5mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-008T-90 (Titanium)

F3-008-45 Castroviejo Suturing Forceps 0.5mm Angled

- Angled jaws, 1 x 2 Teeth
- 0.5mm, with 6mm platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-008T-45 (Titanium)

F3-017-90 Castroviejo Colibri Forceps 0.12mm

- Colibri style shafts, 90° tip
- 0.12mm with 6mm platform, 1 x 2 Teeth
- Overall length: 74mm, 2.9 inches
- Stainless Steel
- F3-017T-90 (Titanium)



F3-017-45 Castroviejo Colibri Forceps 0.12mm

- Colibri style shafts, 45° tip
- 0.12mm with 6mm platform, 1 x 2 Teeth
- Overall length: 74mm, 2.9 inches
- Stainless Steel
- F3-017T-45 (Titanium)



F3-403 Pierse Notched Forceps

- Straight shafts
- 0.3mm hole, Pierse tips with 5mm platform
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-403T (Titanium)



F3-551T St. Martin Forceps Straight

- Straight shafts, 1 x 2 Teeth
- 5mm tying platforms and 0.3mm wide
- Overall length: 87mm, 3.4 inches
- Titanium
- F3-551 (Stainless Steel)



F3-553 St. Martin Forceps Curved

- Curved shafts, 1 x 2 Teeth
- 5mm tying platforms and 0.3mm Tip
- Overall length: 87mm, 3.4 inches
- Stainless Steel
- F3-553T (Titanium)



F3-677 Polack Double Corneal Forceps Colibri Style

- Two 2.75mm long tips serrated 1mm with 1 x 2 teeth
- Serrated handle with dull finish
- Overall length: 72mm, 2.8 inches
- Stainless Steel



F3-715 Suturing Bishop-Harmon Tissue Forceps Delicate Straight

- Straight shafts with 1x 2 teeth, 0.5mm Tip
- Flat handle with holes and polished finish
- Delicate teeth with 3mm platform
- Overall length: 86mm, 3.4 inches
- Stainless Steel



F3-716 Suturing Bishop-Harmon Tissue Forceps Fine Curved

- Curved shafts with 1 x 2 teeth
- With 3mm platform, 0.5mm Tip
- Flat handle with holes and polished finish
- Overall length: 85mm, 3.3 inches
- Stainless Steel



F3-130 Capsule Gimbel Capsulorhexis Forceps 3.0mm Incision 10mm Shafts

- 10mm straight shafts with 6mm platform
- Light tension
- Design to fit through 3.0mm incision
- Flat handle with dull finish
- Overall length: 80mm, 3.2inches
- Stainless Steel



F3-208 Bipolar Castroviejo Curved Capsule Forceps

- 10mm curved shafts
- 2.5mm cups angled 50 degree
- Serrated cross action handle with polished finish
- Overall length: 104mm, 4.1inches
- Stainless Steel



F3-215T Mic Small Incision Forceps 1.8mm Incision 14mm Shafts

- 14mm straight shafts with 4mm platform
- Design to fit through 1.8mm incision
- Overall length: 85mm, 3.3inches
- Titanium



F3-223T Utrata Capsulorhexis 2.8mm Incision 12mm Shafts

- 12mm straight shafts with 4mm platform
- Design to fit through 2.8mm incision
- Overall length: 125mm, 4.9inches
- Titanium
- F3-223 (Stainless Steel)




F3-207 Utrata Capsulorhexis 12mm Shafts

- 12mm straight shafts with 6mm platform
- Overall length: 105mm, 4.13 inches
- Stainless Steel
- F3-207T (Titanium)




F3-211 Utrata Capsulorhexis 1.5mm Incision

- 12mm curved shafts with platform
- Maximum jaws opening at 1.75mm
- Design to fit through 1.5mm incision
- Overall length: 110mm, 4.3inches
- Stainless Steel
- F3-211T (Titanium)



F3-217 Utrata Capsulorhexis 12mm Shafts Curved

- 12mm curved shafts with platform
- Overall length: 110mm, 4.3inches
- Stainless Steel
- F3-217T (Titanium)



F3-227T Micro Incision Forceps 1.8mm Incision 12mm Shafts Straight

- 12mm straight shafts with 6mm platform
- Design to fit through 1.8mm incision
- Small blunt tips insures easy entry
- Overall length: 113mm, 4.4 inches
- Titanium
- F3-227 (Stainless Steel)



F3-251T Mic Small Incision Forceps 1.0mm Incision

- Tapered shafts
- Design to fit through a 1mm incision, 20G
- Overall length: 132mm, 5.6inches
- Titanium
- F3-253T-23 (23G)
- F3-251T-25 (25G)
- F3-251T-27 (27G)



F3-260 Capsule Lehner-Utrata Capsulorhexis Forceps

- 10mm curved shafts
- Utrata-like fine angled tips
- Flat handle with dull finish and platform
- Overall length: 80mm, 3.2inches
- Stainless Steel



F3-288T Micro Incision Forceps 1.8mm Incision 12mm Shafts Curved

- 12mm curved shafts
- Design to fit through 1.8mm incision
- Small blunt tips insures easy entry
- Overall length: 113mm, 4.4 inches
- Titanium
- F3-288 (Stainless Steel)




F3-300T Inamura Capsulorhexis Forceps

- 11mm curved shafts, clear marking on tips
- Cross action flat handle
- Design to fit through 1.8mm incision
- Small blunt tips insures easy entry
- Overall length: 115mm, 4.5 inches
- Titanium




F3-302T Inamura Capsulorhexis Forceps Cross Action

- 10mm curved shafts, clean marking on tips
- Cross action round handle
- Design to fit through 1.5mm incision
- Small blunt tips insures easy entry
- Overall length: 115mm, 4.5 inches
- Titanium



F3-416T Utrata Capsulorhexis 12mm Shafts Straight Angled

- 12mm straight angled shafts
- With 4mm platform
- Overall length: 110mm, 4.3inches
- Titanium
- F3-416 (Stainless Steel)



F3-418 Lehner-Utrata Capsulorhexis Forceps Round Handle

- Utrata-like fine angled tips with gently
- Round handle with dull finish and 4mm platform
- Overall length: 112mm, 4.4inches
- Stainless Steel




F12-002T Lenticule Remove Forceps

- Serrated Tips, 23G
- Curved
- Overall length: 129mm, 5 inches
- Titanium




F3-105T McPherson Tying Forceps 5mm Straight

- Straight shafts
- With 5mm tying platform
- Overall length: 82mm, 3.2 inches
- Titanium
- F3-105 (Stainless Steel)




F3-109T McPherson Tying Forceps 8mm Angled

- Angled shafts
- With 8mm tying platform
- Overall length: 105mm, 4.1 inches
- Titanium
- F3-109 (Stainless Steel)




F3-111 Bonn Tying Forceps 6mm Straight

- Straight shafts
- With 6mm tying platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-111T (Titanium)




F3-112 Bonn Tying Forceps 6mm Curved

- Curved shafts
- With 6mm tying platform
- Overall length: 105mm, 4.1 inches
- Stainless Steel
- F3-112T (Titanium)




F3-632 McPherson Tying Forceps 5mm Angled Short Handle

- Angled shafts
- With 5mm tying platform
- Overall length: 81mm, 3.2 inches
- Stainless Steel
- F3-632T (Titanium)




F3-634 Tying Shepard Forceps 5mm Curved

- Curved shafts
- With 5mm tying platform
- Wide serrated handle
- Overall length: 108mm, 4.3 inches
- Stainless Steel
- F3-634T (Titanium)




F3-636 Tying Shepard Forceps 5mm Angled

- Straight Shafts
- With 5mm tying platform
- Wide serrated handle
- Overall length: 107mm, 4.2 inches
- Stainless Steel
- F3-636T (Titanium)



F3-650 McPherson Tying Forceps 6mm Straight

- Straight shafts
- With 6mm tying platform
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-650T (Titanium)



F3-655 Harms-Tubingen Tying Forceps 6mm Curved

- Curved shafts
- With 6mm tying platform
- Serrated handle with polished finish
- Overall length: 101mm, 4.0 inches
- Stainless Steel
- F3-655T (Titanium)



F3-656 Tying Forceps 7mm Straight

- Straight Shafts
- With 7mm tying platform
- Flat handle
- Overall length: 102mm, 4.0 inches
- Stainless Steel
- F3-656T (Titanium)



F3-670 Tennant Tying Forceps Straight Shafts with Platform

- Straight Shafts
- With 5mm tying platform
- Round knurled handle
- Overall length: 102mm, 4.0 inches
- Stainless Steel
- F3-670T (Titanium)




F3-675T Tennant Tying Forceps Curved Shafts with Platform

- Curved Shafts
- With 5mm tying platform
- Round knurled handle
- Overall length: 103mm, 4.0 inches
- Titanium
- F3-675 (Stainless Steel)



F3-676T Tissue Machat Forceps

- 9mm angled shafts
- With 5mm tying platform
- Smooth tip forceps and smooth the corneal flap in LASIK procedures
- Overall length: 100mm, 3.9 inches
- Titanium





F3-831 Retinal Forceps Diamonized End Gripping

- F3-831-20G (20 gauge)
- F3-831-23G (23 gauge)



F3-835 Retinal Forceps Pic

- F3-835-20G (20 gauge)
- F3-835-23G (23 gauge)
- F3-835-25G (25 gauge)
- F3-835-27G (27 gauge)



F3-837 Retinal Forceps Foreign Body with Cup Jaws

- F3-837-20G (20 gauge)
- F3-837-23G (23 gauge)



F3-838 Retinal Forceps Spring Gripping

- F3-838-20G (20 gauge)
- F3-838-23G (23 gauge)



F3-842 Retinal Forceps ILM

- F3-842-20G (20 gauge)
- F3-842-23G (23 gauge)
- F3-842-25G (25 gauge)
- F3-842-27G (27 gauge)



F3-845 Retinal Forceps Foreign Body

- F3-845-20G (20 gauge)
- F3-845-23G (23 gauge)



F3-849 Retinal Forceps Subretinal

- F3-849-20G (20 gauge)
- F3-849-23G (23 gauge)
- F3-849-25G (25 gauge)
- F3-849-27G (27 gauge)



Tissue • Utility • Iris • Fixation • Prechopper Forceps

F3-233T Iris Jeweler Type Forceps Fine Tying Platform

- Straight shafts with pointed tips
- 6mm tying platform
- Fine pointed tips
- Overall length: 120mm, 4.7 inches
- Titanium



F3-234 Iris Jeweler Type Forceps Fine Curved

- Curved shafts with fine pointed tips
- Smooth handle with dull finish
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F3-236 Bipolar Jeweler

- Straight shafts with micro tips
- Dull finish
- Smooth handle with dull finish
- Overall length: 119mm, 4.7 inches
- Stainless Steel



F3-276 Leduc Meibomian Forceps

- Gentle paddle shaped tips
- Overall length: 90mm, 3.5 inches
- Stainless Steel



F3-277 Moorfields Forceps

- Serrated jaws
- 1mm x 1.25mm tips
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F3-521T Iris Forceps Straight

- Straight shafts with serrated tips
- Overall length: 106mm, 4.2 inches
- Titanium
- F3-521 (Stainless Steel)



F3-523 Iris Forceps Angled

- Angled shafts with serrated tips
- Overall length: 106mm, 4.2 inches
- Stainless Steel
- F3-523T (Titanium)



F3-543 Akahoshi Universal Prechopper Angled


- 22mm angled shafts
- Sharp point tips
- Overall length: 120mm, 4.7 inches
- Stainless Steel



Tissue • Utility • Iris • Fixation • Prechopper Forceps



F3-545 Akahoshi Universal Prechopper Straight

- Straight shafts with sharp point tips
- 22mm from tip to bend
- Overall length: 120mm, 4.7 inches
- Stainless Steel




F3-549 Akahoshi Combo Prechopper

- 23mm straight shafts with crescent tips
- Overall length: 120mm, 4.7 inches
- Stainless Steel
- F3-549T (Titanium)




F3-550 Femto Prechopper

- Overall length: 120mm, 4.7 inches
- Stainless Steel



F3-566 Cilia Barraquer Forceps

- Slightly curved with 6mm platform
- Serrated handle with polished finish
- Overall length: 109mm, 4.3 inches
- Stainless Steel
- F3-566T (Titanium)




F3-701 Dressing Forceps Straight with Serrations Delicate

- Straight shafts
- Overall length: 100mm, 3.9 inches
- Stainless Steel




F3-702 Fixation Lester Fixation Forceps 2 x 3

- Straight shafts and serrated handle with polished finish
- Fine 2 x 3 teeth
- Overall length: 96mm, 3.8 inches
- Stainless Steel



F3-703 Dressing Forceps Curved with Serrations Delicate

- Curved shafts with 0.5mm distance serrated tips
- Overall length: 100mm, 3.9 inches
- Stainless Steel



F3-704 Fixation Lester Moody Curved Fixation Forceps with Lock

- Curved shafts with 0.5mm, 1 x 2 teeth
- Serrated handle with thumb catch lock and polished finish
- Overall length: 108mm, 4.3 inches
- Stainless Steel



Tissue • Utility • Iris • Fixation • Prechopper Forceps

F3-741T Conjunctiva Forceps

- Straight handle with hole on tip
- With 6mm tying platform
- Overall length: 82mm, 3.2 inches
- Titanium
- F3-741(Stainless Steel)



F3-761 Fixation Graefe Fixation Forceps Narrow 3.5mm

- Straight jaws with fine teeth
- Serrated handle with polished finish
- Thumb without catch lock
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F3-765 Fixation Graefe Fixation Forceps Narrow 3.5mm with Lock

- Straight jaws with fine teeth
- Serrated handle with polished finish
- Thumb catch lock
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F3-771 Chalazion Forceps 20mm x 13mm

- Chalazion forceps oval
- Large, 20mm x 13mm inside dimensions
- Overall length: 90mm, 3.5 inches
- Stainless Steel



F3-773 Chalazion Forceps 14mm x 11mm

- Chalazion forceps oval
- Medium, 14mm x 11mm inside dimensions
- Overall length: 90mm, 3.5 inches
- Stainless Steel



F3-775 Chalazion Forceps 8mm

- Chalazion forceps oval
- Small, Round, 8mm x 8mm inside dimensions
- Overall length: 90mm, 3.5 inches
- Stainless Steel



F3-801 Serrefine

- Serrated jaws
- Small, 11 / 20 (38mm) long
- Overall length: 70mm, 2.8 inches
- Stainless Steel



Lens Implantation • Lasik • Lasek • Flap Forceps

F5-057 Halstead Curved Hemostatic Mosquito Forceps

- Curved 18mm serrated jaws
- Right handle with ratchet lock
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F5-055 (Straight)



F3-163 Donnor Membrane Holder / Dmek

- Adapt for DMEK surgery with broad grasping
- the endothelium membrane of the donor
- Overall length: 100mm, 3.9 inches
- Stainless Steel



F3-301 Intraocular Lens Forceps Angled

- 10mm angled shafts
- With 1mm diameter jaws
- Overall length: 100mm, 3.9 inches
- Stainless Steel
- F3-301T (Titanium)



F3-303 Intraocular Lens Forceps Curved

- Curved shafts with 1mm diameter jaws
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F3-303T (Titanium)



F3-305T Intraocular Lens Forceps Angled

- 9mm angled tips
- Overall length: 125mm, 4.9 inches
- Titanium
- F3-305 (Stainless Steel)



F3-306 Intraocular Lens Implant Forceps

- Provides grasp of the lens optic or segment
- Thin instrument shaft
- Overall length: 120mm, 4.8 inches
- Stainless Steel



F3-308 Intraocular Lens IOL Forceps

- Used with pre-positioned,
- Winged cartridge IOL delivery systems
- Work well with both acrylic and silicone implants
- Overall length: 94mm, 3.7 inches
- Stainless Steel



F3-321T Lens Implantation Forceps

- 9mm angled shafts
- 1mm x 2mm dentation jaws
- Overall length: 105mm, 4.1 inches
- Titanium



Lens Implantation • Lasik • Lasek • Flap Forceps

F3-323T Intraocular Lens Universal Folding Forceps

- Idea for lens diameters from 5.5mm and greater
- Tip is for tucking the haptic between the lens halves
- Overall length: 105mm, 4.1 inches
- Titanium



F3-325T Lens Forceps with Lock

- 20mm angled shafts
- With lock facilitates
- Design to fit through 2mm incision
- Overall length: 110mm, 4.3 inches
- Titanium




F3-327T Lens Implantation Forceps

- 20mm angled shafts
- Design to fit through 2mm incision
- Overall length: 110mm, 4.3 inches
- Titanium



F3-329T Lens Buratto Implantation Forceps Angled

- Angled shafts
- For HYDROVIEW lens
- Overall length: 105mm, 4.1 inches
- Titanium



F3-330 Intraocular Lens Shepard Curved Forceps

- Curved shafts
- Serrated 1.5mm upper and 2mm lower jaws
- Handle with slide lock and polished finish
- Overall length: 91mm, 3.6 inches
- Stainless Steel



F3-332T Intraocular Shepard IOL Insertion Forceps

- 18mm curved shafts with 1.91mm wide grasping tip
- Flat, serrated cross action handle with dull finish
- Overall length: 108mm, 4.3 inches
- Titanium



F3-501 Strabismus Forceps 4 x 4 Teeth

- 10mm angled shafts, 90 degree
- 1mm teeth on one jaw fit into holes on the other jaw
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-501-L (Left, 4 x 4 Teeth)
- F3-501-R (Right, 4 x 4 Teeth)



F3-555 Perone Lasik Flap Forceps

- Colibri style tips with dia. 1.7mm ring
- 0.3mm Hole
- Overall length: 105mm, 4.1 inches
- Stainless Steel



Lens Implantation • Lasik • Lasek • Flap Forceps

F3-559 Retinal Sclera Plug Forceps

- 10mm curved shafts with 1.5mm wide jaws
- Overall length: 115mm, 4.5 inches
- Stainless Steel
- F3-559T (Titanium)



F3-661T Lens Forceps

- Shafts angled 45 degree
- Curved grooved bladed for holding the landers vitrectomy lenses
- Landers vitrectomy
- Overall length: 115mm, 4.5 inches
- Titanium



F3-688T Artisan Lens Implantation Forceps

- 7mm angled shafts
- Overall length: 124mm, 4.9 inches
- Titanium



F3-788T Trocar Ports Grasping Forceps

- 23 gauge
- For Easy Removal of Trocar Port From Sclera
- Overall length: 113mm, 4.4 inches
- Titanium



F3-802T-R Artiflex Lens Implantation Forceps Right

- Overall length: 107mm, 4.2 inches
- Titanium



F3-802T-L Artiflex Lens Implantation Forceps Left

- Overall length: 107mm, 4.2 inches
- Titanium



F3-505 Muscle Clamp Curved Shafts

- Curved shafts
- 0.25mm teeth without tying platform
- Wide serrated handle with thumb catch lock
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-505-L (Left)
- F3-505-R (Right)



F3-509 Muscle Clamp Straight Shafts

- Straight shafts
- 0.25mm teeth without tying platform
- Wide serrated handle with thumb catch lock
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F3-513 Jameson Muscle Forceps 4 x 4 Teeth

- 10mm angled shafts, 90 degree
- Teeth on one jaw fit into holes on the other jaws
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-513-L (Left, 4 x 4 Teeth)
- F3-513-R (Right, 4 x 4 Teeth)



F3-514 Jameson Muscle Forceps Child size

- 9mm angled shafts, 45 degree
- 1mm teeth fit into the other jaw
- Slide lock and serrated handle with polished finish
- Overall length: 99mm, 3.9 inches
- Stainless Steel
- F3-514-L (Left)
- F3-514-R (Right)



F3-517 Jameson Muscle Forceps 3 x 3 Teeth

- 7mm angled shafts, 90 degree
- Teeth on one jaw fit into holes on the other jaws
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F3-517-L (Left, 3 x 3 Teeth)
- F3-517-R (Right, 3 x 3 Teeth)



F3-588T Muscle Clamp

- Putterman clamp
- T mode
- Overall length: 94mm, 3.7 inches
- Titanium
- F3-588 (Stainless Steel)



Scissors

Retinal

Corneal

Vannas


Iris


Tenotomy

Eye

Lens

Dalk





F3-751 Retinal Scissors Straight Blunt Tip

- F3-751-20G (20 gauge)
- F3-751-23G (23 gauge)



F3-752 Retinal Scissors Vertical 70 degree

- F3-752-20G (20 gauge)
- F3-752-23G (23 gauge)



F3-756 Retinal Scissors Horizontal 55 degree

- F3-756-20G (20 gauge)
- F3-756-23G (23 gauge)



F3-758 Retinal Scissors Angled 45 degree

- F3-758-20G (20 gauge)
- F3-758-23G (23 gauge)



F3-759 Retinal Scissors Curved Subretinal

- F3-759-20G (20 gauge)
- F3-759-23G (23 gauge)



Blunt Tip



Sharp Tip



F4-001 Universal Corneal Scissors 14mm Blunt Tip Curved

- Medium blades with blunt tips
- 9mm cut length, 14mm from mid-screw to tip
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F4-001T (Titanium)

F4-003 Universal Corneal Scissors 14mm Sharp Tip Curved

- Medium blades with sharp tips
- 9mm cut length, 14mm from mid-screw to tip
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F4-003T (Titanium)

F4-005 Universal Corneal Scissors 16mm Blunt Tip Curved

- Medium blades with blunt tips
- 10.5mm cut length, 16mm from mid-screw to tip
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F4-005T (Titanium)

F4-015 Universal Corneal Scissors 12mm Sharp Tip Curved

- Medium blades with sharp tips
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F4-015T (Titanium)

F4-007 Universal Corneal Scissors 14mm Blunt Tip Curved

- Medium blades with blunt tips
- 9mm cut length, 14mm from mid-screw to tip
- Overall length: 108mm, 4.3 inches
- Stainless Steel
- F4-007T (Titanium)

F4-011 Universal Corneal Scissors 12mm Blunt Tip Curved

- Medium blades with blunt tips
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F4-011T (Titanium)

F4-017 Corneal Scissors

- Corneoscleral scissors
- Small blades
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 108mm, 4.3 inches
- Stainless Steel
- F4-017-L (Left curved)
- F4-017-R (Right curved)



F4-055 Castroviejo Curved Corneal Scissors

- Curved, blunt tips with flat handle
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F4-056 Castroviejo Curved Corneal Scissors

- Curved, sharp pointed tips with round handle
- 10.5mm cut length, 16mm from mid-screw to tip
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F4-101 Corneal Scissors

- Universal miniature model with curved blunt tips
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 100mm, 4.0 inches
- Stainless Steel
- F4-101T (Titanium)





Straight



Curved



F4-031 Vannas Scissors 9mm Straight

- Straight and sharp pointed tips
- 6mm cut length, 9mm from mid-screw to tip
- Overall length: 85mm, 3.4 inches
- Stainless Steel
- F4-031T (Titanium)

F4-033 Vannas Scissors 9mm Curved

- Curved and sharp pointed tips
- 6mm cut length, 9mm from mid-screw to tip
- Overall length: 85mm, 3.4 inches
- Stainless Steel
- F4-033T (Titanium)

F4-035 Vannas Scissors 7mm Straight

- Straight and sharp pointed tips
- 7mm cut length, 12mm from mid-screw to tip
- Overall length: 105mm, 4.13 inches
- Stainless Steel
- F4-035T (Titanium)

F4-037 Vannas Scissors 7mm Curved

- Curved and sharp pointed tips
- 7mm cut length, 12mm from mid-screw to tip
- Overall length: 105mm, 4.13 inches
- Stainless Steel
- F4-037T (Titanium)

F4-041 Vannas Capsulotomy Scissors 12mm Straight

- Straight and sharp pointed tips
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F4-041T (Titanium)

F4-043 Vannas Capsulotomy Scissors 12mm Curved

- Curved and sharp pointed tips
- 8mm cut length, 12mm from mid-screw to tip
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F4-043T (Titanium)

F4-045 Vannas Capsulotomy Scissors 16mm Straight

- Straight and sharp pointed tips
- 10.5mm cut length, 16mm from mid-screw to tip
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F4-045T (Titanium)

F4-047 Vannas Capsulotomy Scissors 16mm Curved

- Curved and sharp pointed tips
- 10.5mm cut length, 16mm from mid-screw to tip
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F4-047T (Titanium)

F4-051 Tenotomy Scissors 20mm Sharp Tip Curved

- Punctal scissors
- One conically shaped blade with micro grooves and textured surface
- 13.5mm cut length, from mid-screw to tip 20mm
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F4-053 Tenotomy Scissors 20mm Blunt Tip Curved

- Blunt tips
- Medium blades
- 13.5mm cut length, from mid-screw to tip 20mm
- Overall length: 120mm, 4.7 inches
- Stainless Steel



F4-057 Bonn Model Eye Scissors

- Dull straight
- Pointed tips
- 13.5mm cut length, from mid-screw to tip 20mm
- Overall length: 85mm, 3.4 inches
- Stainless Steel
- F4-059 (Pointed tips curved)



F4-072 Eye Scissors Straight

- Pointed tips
- Standard
- 22mm cut length, from mid-screw to tip 33mm
- Overall length: 100mm, 4.0 inches
- Stainless Steel
- F4-074 (Blunt tips straight)



F4-076 Eye Scissors Curved

- Curved tips
- Standard
- 22mm cut length, from mid-screw to tip 33mm
- Overall length: 100mm, 4.0 inches
- Stainless Steel
- F4-078 (Blunt tips curved)



F4-080 Eye Scissors Straight

- Pointed tips
- Standard
- 20mm cut length, from mid-screw to tip 30mm
- Overall length: 138mm, 5.4 inches
- Stainless Steel



F4-081 Enucleation Scissors Medium Light Curved

- Medium light curved
- Blunt tips
- 35mm cut length, from mid-screw to tip 52mm
- Overall length: 125mm, 5.0 inches
- Stainless Steel
- F4-083 (Large light curved)



F4-111 Barraquer Iris Scissors

- Cross action flat serrated handle
- Delicate angled 9mm blades, sharp pointed tips
- 12mm cut length, from mid-screw to tip 17mm
- Overall length: 85mm, 3.4 inches
- Titanium
- F4-111T (Titanium)



F4-085 Stevens Tenotomy Scissors Straight

- Straight blades and blunt tips
- 17mm cut length, from mid-screw to tip 30mm
- Polished finish on blades
- Overall length: 115mm, 4.5 inches
- Stainless Steel



F4-086 Stevens Tenotomy Scissors Curved

- Curved blades and blunt tips
- 17mm cut length, from mid-screw to tip 30mm
- Polished finish on blades
- Overall length: 115mm, 4.5 inches
- Stainless Steel



F4-091T Lens Cutter

- Designed to cut and remove all types of foldable intraocular lens through a small incision
- 8mm cut length, from mid-screw to tip 12mm
- Overall length: 125mm, 5.0 inches
- Titanium



F4-361T Dalk Scissors

- Used to remove the superficial layers of corneal stroma.
- 7mm cut length, from mid-screw to tip 10.5mm
- Overall length: 106mm, 4.2 inches
- Titanium
- F4-361 (Stainless Steel)



Needle Holder

Needle Holders



F5-001 Barraquer Needle Holder 10mm Round Handle

- Fine 10mm gently curved jaws
- Round knurled handle
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F5-001T (Titanium)



F5-003T Barraquer Needle Holder 10mm Jaws Curved Round Handle

- Delicate 10mm gently curved jaws
- Cross action round knurled handle
- Overall length: 125mm, 4.9 inches
- Titanium
- F5-003 (Stainless Steel)



F5-005 Barraquer Needle Holder 10mm Jaws Curved

- Delicate 10mm gently curved jaws
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F5-005T (Titanium)



F5-007T Barraquer Needle Holder with Lock Round Handle

- Delicate 10mm gently curved jaws
- Round knurled handle with solid lock
- Overall length: 125mm, 4.9 inches
- Titanium
- F5-007 (Stainless Steel)



F5-009 Barraquer Needle Holder 8mm Strong Jaws Curved Round Handle

- Delicate 8mm gently curved jaws
- Cross action round knurled handle
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F5-009T (Titanium)



F5-011T Barraquer Needle Holder 9mm Fine Jaws Round Handle

- Fine 9mm gently curved jaws
- Round knurled handle
- Overall length: 100mm, 3.9 inches
- Titanium
- F5-011 (Stainless Steel)



F5-013 Barraquer Needle Holder 10mm Jaws Curved Round Handle

- Delicate 10mm gently curved jaws
- Round knurled handle
- Overall length: 125mm, 4.9 inches
- Stainless Steel
- F5-013T (Titanium)



F5-015T Castroviejo Needle Holder Strong Jaws Flat Handle

- Delicate 10mm gently curved jaws
- Flat handle
- Overall length: 125mm, 4.9 inches
- Titanium
- F5-015 (Stainless Steel)



F5-017 Castroviejo Needle Holder with Lock Flat Handle

- 10mm standard curved jaws
- Flat handle with lock
- Overall length: 135mm, 5.3 inches
- Stainless Steel
- F5-017T (Titanium)



F5-019 Castroviejo Needle Holder 15mm Flat Handle

- 15mm standard curved jaws
- Flat handle
- Overall length: 140mm, 5.5 inches
- Stainless Steel
- F5-019T (Titanium)



F5-021 Castroviejo Needle Holder 10mm Delicate Jaws Flat Handle

- Delicate 10mm gently curved jaws
- Cross action flat handle
- Overall length: 135mm, 5.3 inches
- Stainless Steel
- F5-021T (Titanium)



F5-023 Needle Holder Scissors 16mm Fine Jaws Curved Round Handle

- Fine 16mm gently curved jaws
- With 9mm blades from mid-screw to jaws
- Overall length: 127mm, 5.0 inches
- Stainless Steel
- F5-025 (Straight Fine Jaws)



F5-031 Kalt Needle Holder

- Medium 10mm finely serrated straight jaws
- Overall length: 130mm, 5.1 inches
- Stainless Steel



F5-061 Baby Jones Towel Clamp

- Medium 23mm finely serrated angled jaws
- Overall length: 90mm, 3.5 inches
- Stainless Steel



Caliper and Marker

Caliper
Marker



F6-001 Fixation Ring

- 16mm smooth bulbous tip
- Overall length: 40mm, 1.6 inches
- Stainless Steel
- F6-003 (18mm tips, Overall length: 43mm)
- F6-005 (20mm tips, Overall length: 45mm)
- F6-007 (22mm tips, Overall length: 48mm)



F6-011T Castroviejo Caliper Straight Shafts

- Measure from 0mm to 20mm in 1mm increment
- Scale reads on one side
- Straight
- Overall length: 80mm, 3.2 inches
- Titanium
- F6-011 (Stainless Steel)



F6-013 Castroviejo Caliper Curved Shafts

- Measure from 0mm to 20mm in 1mm increment
- Scale reads on one side
- Curved 45 degree shafts
- Overall length: 80mm, 3.2 inches
- Stainless Steel
- F6-013T (Titanium)




F6-021T Fixed Caliper 3.0mm / 3.5mm

- 3.0mm / 3.5mm
- Fine tips mark sclerotomy sites for vitreous surgery
- Dull finish
- Overall length: 70mm, 2.8 inches
- Titanium
- F6-021 (Stainless Steel)



F6-023 Fixed Caliper 3.5mm / 4.0mm

- 3.5mm / 4.0mm
- Fine tips mark sclerotomy sites for vitreous surgery
- Dull finish
- Overall length: 70mm, 2.8 inches
- Stainless Steel
- F6-023T (Titanium)



F6-030T Axis Marker

- Two line marker
- Used with Mendez Degree Gauge
- Overall length: 81mm, 3.4 inches
- Titanium



F6-035T Low Profile Corneal Marker 8 Line

- 8 Line
- Round serrated handle
- Osher-neuman low profile
- Overall length: 136mm, 5.4 inches
- Titanium
- F6-031T (4 Line)
- F6-033T (6 Line)



F6-036 Toric Axis Marker

- Overall length: 100mm, 3.9 inches
- Stainless Steel



F6-037 Toric Reference Ring Marker

- Overall length: 105mm, 4.1 inches
- Stainless Steel



F6-039T Mendez Marker

- Calibrated every 10 degrees
- From 0 degree to 180 degree
- Internal diameter 13mm, External diameter 19mm
- Overall length: 136mm, 5.4 inches
- Titanium



F6-040 McDonald Marker 4 Line

- Used in mini-RK to create four incision
- Inner ring measures in 4mm diameter, outer ring measures 5mm in diameter
- Marker can be used in 4, 4.5, and 5 optical zone
- Overall length: 107mm, 4.2 inches
- Stainless Steel



F6-042 Osher-Neumann Corneal Marker 8 Line

- 8 Line
- Marks radial lines is 4mm, Inner diameter of blades is 4.5mm
- Overall length: 109mm, 4.3 inches
- Stainless Steel
- F6-038 (12 Line)
- F6-046 (16 Line)



F6-050 Toric Bubble Marker

- The marker is pressed against the corneal surface to provide marks at the horizontal axis and at 90
- Overall length: 99mm, 3.9 inches
- Stainless Steel



F6-052 Wallace Capsulorhexis Gauge 6mm

- 6mm diameter
- To gauge the diameter of the anterior capsular opening
- Overall length: 81mm, 3.2 inches
- Stainless Steel



F6-054T Whitman LRI Marker

- Astigmatic dial marked in degrees
- Provides two large radial blades
- Overall length: 81mm, 3.2 inches
- Titanium
- F6-054 (Stainless Steel)



F6-056 Whitman Double- Ended LRI Marker

- The two astigmatic dials cover a wide range of dioptric treatments
- Overall length: 122mm, 4.8 inches
- Stainless Steel



F6-770T Corneal Marker 8.5mm / 9mm

- One end 8.5mm diameter, One end 9mm diameter
- Provide a template for scoring and stripping during the DSEK procedure
- Overall length: 132mm, 5.2 inches
- Stainless Steel



F6-780 Gulani LASIK Marker

- 3.5mm and 4mm intersecting circles to mark corneal flap for LASIK procedures
- Round handle
- Overall length: 119mm, 4.7 inches
- Stainless Steel



F6-912 MacRae mini LASEK / PRK Well 6.5mm

- 6.5mm diameter
- Used for PRK or LASEK
- Remove less epithelium
- Overall length: 81mm, 3.2 inches
- Stainless Steel
- F6-914 (7mm diameter)



F9-715 Fine Thornton Fixation Ring

- Globe fixation ring with swivel handle 16mm diameter ring multiple blunt tip
- With swivel
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F9-717 Fine Thornton Fixation Ring

- 16mm internal diameter rotatable ring with atraumatic teeth provides excellent stabilization and adequate pressure on globe
- With swivel
- Overall length: 110mm, 4.3 inches
- Stainless Steel



F9-751 Optic Zone Marker

- 3mm diameter
- Height of ring 2.3mm
- Overall length: 122mm, 4.8 inches
- Stainless Steel
- F9-753 (4mm diameter)
- F9-755 (5mm diameter)
- F9-757 (6mm diameter)
- F9-759 (7mm diameter)
- F9-761 (8mm diameter)



F9-849 Jaeger Lid Plate

- 18mm-22mm curved wide blades
- Overall length: 100mm, 3.9 inches
- Stainless Steel



F9-900T Lasek Trephine

- Sharp ring in 0.5mm increments
- 7.5mm diameter
- Overall length: 121mm, 4.8 inches
- Titanium
- F9-902T (8.0mm diameter)
- F9-904T (8.5mm diameter)
- F9-906T (9.0mm diameter)



F9-901T Lasek Funnel

- Blund ring in 0.5mm increments
- 7.5mm diameter
- Overall length: 123mm, 4.8 inches
- Titanium
- F9-903T (8.0mm diameter)
- F9-905T (8.5mm diameter)
- F9-907T (9.0mm diameter)



F12-001 Fixation Ring

- Diameter 13mm
- Overall length: 107mm, 4.2 inches
- Stainless Steel



07


Focus[®]
INSTRUMENTS


Ronguer Punch Trephine

Rongeur

Scleral Punch

Trephine





F7-001 Kerrison Rongeur 2mm x 3mm 90 degree

- 2mm x 3mm 90 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel

F7-003 Kerrison Rongeur 3mm x 4.5mm 90 degree

- 3mm x 4.5mm 90 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel

F7-005 Kerrison Rongeur 4mm x 3.5mm 90 degree

- 4mm x 3.5mm 90 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel

F7-007 Kerrison Rongeur 2mm x 3mm 100 degree


- 2mm x 3mm 100 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel

F7-009 Kerrison Rongeur 3mm x 4.5mm 100 degree

- 3mm x 4.5mm 100 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel

F7-011 Kerrison Rongeur 4mm x 3.5mm 100 degree

- 4mm x 3.5mm 100 degree
- Sharp blade tips for lacrimal surgery
- Polished finish
- Stainless Steel



Straight



F7-021 Descemet Membrane Punch 1.0mm Straight

- Straight 1.0mm
- 0.9mm x 0.45mm
- Designed to produce a clean 1.0mm round hole without tissue tags
- Smooth squeeze action handle with holds
- Overall length: 125mm, 4.9 inches
- Stainless Steel

F7-023 Descemet Membrane Punch 1.5mm Straight

- Straight 1.5mm
- 2mm x 1mm
- Designed to produce a clean 1.5mm round hole without tissue tags
- Smooth squeeze action handle with holds
- Overall length: 125mm, 4.9 inches
- Stainless Steel

F7-025 Descemet Membrane Punch 1.0mm Curved

- Curved 1.0mm
- 0.9mm x 0.45mm
- Designed to produce a clean 1.0mm round hole without tissue tags
- Smooth squeeze action handle with holds
- Overall length: 125mm, 4.9 inches
- Stainless Steel

F7-027 Descemet Membrane Punch 1.5mm Curved

- Curved 1.5mm
- 2mm x 1mm
- Designed to produce a clean 1.5mm round hole without tissue tags
- Smooth squeeze action handle with holds
- Overall length: 125mm, 4.9 inches
- Stainless Steel



F7-031 Trephine 5mm

- 5mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-033 Trephine 5.25mm

- 5.25mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-035 Trephine 5.55mm

- 5.55mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-037 Trephine 5.75mm

- 5.75mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-039 Trephine 6mm

- 6mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-041 Trephine 6.25mm

- 6.25mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-043 Trephine 6.5mm

- 6.5mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-045 Trephine 6.75mm

- 6.75mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-047 Trephine 7mm

- 7mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-049 Trephine 7.25mm

- 7.25mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-051 Trephine 7.5mm

- 7.5mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-053 Trephine 7.75mm

- 7.75mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-055 Trephine 8mm

- 8mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-057 Trephine 8.25mm

- 8.25mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-059 Trephine 8.55mm

- 8.55mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-061 Trephine 8.75mm

- 8.75mm in 0.25mm increments
- Dull finish
- Stainless Steel

F7-063 Trephine 9mm


- 9mm in 0.25mm increments
- Dull finish
- Stainless Steel

08

Focus[®]
INSTRUMENTS

Hooks Lacrimal Retractor

Hooks
Lacrimal
Retractor





F8-061 Jameson Strabismus Hook 6mm

- Muscle hook
- 6mm tip to curve
- Overall length: 123mm, 4.8 inches
- Stainless Steel
- F8-063 (8mm)
- F8-065 (10mm)
- F8-067 (12mm)



F8-069 Jameson Strabismus Hook 6mm with Hole

- Muscle hook
- 6mm tip to curve
- Overall length: 123mm, 4.8 inches
- Stainless Steel
- F8-071 (8mm)
- F8-073 (10mm)
- F8-075 (12mm)



F8-077 Jameson Strabismus Hook 8mm

- Muscle hook
- 8mm tip to curve
- Extra delicate with blunt spatulated tip
- Overall length: 123mm, 4.8 inches
- Stainless Steel
- F8-079 (10mm)
- F8-081 (12mm)



F8-083 Jameson Strabismus Hook 10mm

- Muscle hook
- 10mm tip to curve
- Extra delicate with blunt spatulated tip
- Overall length: 123mm, 4.8 inches
- Stainless Steel
- F8-085 (12mm)



F8-089 Jameson Strabismus Hook 6mm

- Extra delicate, teaser hook
- 6mm tip to curve
- Blunt tip
- Overall length: 123mm, 4.8 inches
- Stainless Steel



F8-091 Jameson Strabismus Hook 2.5mm

- Extra delicate, teaser hook
- 2.5mm tip to curve
- Blunt tip
- Overall length: 123mm, 4.8 inches
- Stainless Steel



F8-093 Jameson Strabismus Hook 1.5mm

- Double arm
- 1.5mm wide
- Overall length: 123mm, 4.8 inches
- Stainless Steel



F8-097 Jameson Strabismus Hook 1.75mm

- 1.75mm wide
- Scleral ruler 15mm long with notches in 5mm increments
- single marking tooth at tip
- For measuring and marking the distance from the limbus to the location of the muscle reinsertion
- Overall length: 123mm, 4.8 inches
- Stainless Steel



F8-095 Jameson Strabismus Hook 0.7mm

- Muscle hook
- 0.7mm
- Overall length: 123mm, 4.8 inches
- Stainless Steel



F8-001 Desmarres Lid Retractor

- 8mm wide retractor blade
- Overall length: 128mm, 5.0 inches
- Stainless Steel
- F8-003 (10mm blade)
- F8-005 (12mm blade)



F8-015 Arruga Globe Retractor S Mode Long

- S2 mode long
- 7.5mm width x 3.6mm
- Polished finish
- Overall length: 135mm, 5.3 inches
- Stainless Steel



F8-017 Arruga Globe Retractor S Mode Short

- S2 mode short
- 3.6mm x 9.5mm
- Round knurled handle
- Overall length: 130mm, 5.1 inches
- Stainless Steel



F8-019 Arruga Globe Retractor Z Mode

- Z mode
- Double ended with 8mm and 10mm wide retractor blade between
- Flat handle
- Overall length: 120mm, 4.7 inches
- Stainless Steel
- F8-019T (Titanium)



F8-023T Arruga Globe Retractor Flat

- 8mm x 25mm wide retractor blade
- Flat handle
- Overall length: 135mm, 5.3 inches
- Titanium
- F8-027T (6mm x 25mm)



F8-101 Lid Retractor 6mm / 7.2mm / 9.2mm

- Conjunctiva retractor
- 6mm wide
- Thin solid blade
- Overall length: 123mm, 4.8 inches
- Stainless Steel
- F8-103 (7.2mm wide)
- F8-105 (9.2mm wide)



F2-099 Infusion Lens Loop Hook

- 3mm wide x 8mm long ultra thin loop
- Round knurled handle
- Overall length: 118mm, 4.7 inches
- Stainless Steel
- F2-101 (5mm x 7mm)
- F2-103 (4mm x 6mm)



F2-105 Infusion Lens Loop Spoon

- 3mm wide x 9mm long blunt curette type cup
- Cataract spoon
- Overall length: 130mm, 5.1 inches
- Stainless Steel



F2-171 Lacrimal Dilators

- Long 24mm taper with blunt tip
- Round knurled handle
- Polished finish
- Overall length: 100mm, 3.9 inches
- Stainless Steel
- F2-173 (Medium)
- F2-175 (Short)



F9-801T Schocket Scleral Depressor Straight

- Depressor double straight
- 2x5 mm cylinder tip and
- 2x6 mm bar tip
- Overall length: 133mm, 5.2 inches
- Titanium



F9-803T Schocket Scleral Depressor Curved

- Depressor double curved
- 2x6 mm cylinder tip and
- 2x7 mm bar tip
- Overall length: 133mm, 5.2 inches
- Titanium



F9-809T Schocket Scleral Depressor

- 6mm angled shaft
- 2mm teardrop shaped tip
- Overall length: 120mm, 4.7 inches
- Titanium



F9-811T Schocket Scleral Depressor

- 6mm wide bar and 0.5mm teardrop shape tips
- Overall length: 133mm, 5.2 inches
- Titanium



F9-843 Meyerhoefer Chalazion Curette

- 3.5mm cup
- Straight 30mm shaft
- Flat handle
- Overall length: 130mm, 5.1 inches
- Stainless Steel
- F9-845 (2.5mm cup)
- F9-847 (1.5mm cup)
- F9-843X (1.0mm cup)



F9-851 Pigtail Probe

- 0.6mm diameter pigtail curved probes with suture eyes just back of blunt ends
- 23g
- Overall length: 110mm, 4.3 inches
- Stainless Steel
- F9-853 (0.7mm diameter, 22g)
- F9-855 (0.8mm diameter, 21g)
- F9-857 (0.9mm diameter, 20g)



09

Focus[®]
INSTRUMENTS

Phaco Chopper Manipulator Spatula

Phaco Chopper
Manipulator
Spatula
Polisher





F9-101 Nagahara Nucleus Chopper L-hand 45 degree Angle 1.3mm

- 45 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.3mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-101T (Titanium)

F9-103 Nagahara Nucleus Chopper L-hand 45 degree Angle 1.5mm

- 45 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.5mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-103T (Titanium)

F9-105 Nagahara Nucleus Chopper L-hand 45 degree Angle 1.75mm

- 45 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.75mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-105T (Titanium)

F9-107 Nagahara Nucleus Chopper L-hand 45 degree Curved

- 45 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.5mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-107T (Titanium)

F9-109 Nagahara Nucleus Chopper L-hand 45 degree Curved 1.8mm

- 45 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.8mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-109T (Titanium)



F9-111 Nagahara Nucleus Chopper L-hand 90 degree Angle

- 90 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.3mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-111T (Titanium)

F9-113 Nagahara Nucleus Chopper L-hand 90 degree Angle


- 90 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.5mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-113T (Titanium)

F9-115 Nagahara Nucleus Chopper L-hand 90 degree Angle

- 90 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.8mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-115T (Titanium)

F9-206 Nagahara Nucleus Chopper L-hand 60 degree Angle

- 60 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.3mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel



F9-202 Nagahara Nucleus Chopper L-hand 60 degree Angle

- 60 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.75mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel

F9-204 Nagahara Nucleus Chopper L-hand 60 degree Angle

- 60 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 1.5mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel

F9-130 Nagahara Nucleus Chopper L-hand 60 degree Angle

- 60 degree Nagahara on one end to chop towards phaco handpiece
- Tip length 2.5mm
- Nagahara style
- Overall length: 120mm, 4.7 inches
- Stainless steel

F9-117 Nucleus Chopper B-shape

- B-shape rotator tip
- 11mm angled shaft
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-117T (Titanium)



F9-119T Jewelewicz Rosen Chopper

- L-shape rotator tip
- 10mm angled shaft with pointed tip
- Tip length 1.3mm
- Overall length: 120mm, 4.7 inches
- Titanium
- F9-119 (Stainless Steel)




F9-210 Kavanagh Nucleus Seibel Chopper

- 0.6mm olive shape tip
- Tip length: 1.75mm
- Overall length: 115mm, 4.6 inches
- Stainless steel



F9-216 Kavanagh Nucleus Seibel Chopper

- Using a prechop forceps for dense(4+) nuclear sclerotic cataract
- 10mm angled shaft, 45 degree
- 0.5mm diameter bulbous tip
- Tip length 1.7mm
- Overall length: 123mm, 4.8 inches
- Stainless steel



F9-226T Wallace Double-Ended Guardian & Quickchopper

- Combines the Wallace Guardian Chopper with a quick chopper tip
- Gently angled blunt chopping tip
- Angled shaft for ease of manipulation
- Overall length: 132mm, 5.2 inches
- Titanium



F9-228T Wallace Guardian Chopper

- Gently angled blunt chopping tip
- Angled shafts for ease of manipulation
- Overall length: 112mm, 4.4inches
- Titanium



F9-001 Sinskey Hook

- 10mm angled shaft 45 degree
- 0.2mm hook
- Overall length: 120mm, 4.7 inches
- Stainless Steel
- F8-001T (Titanium)



F9-003 Fenzl Hook

- Fenzl hook
- Overall length: 120mm, 4.7 inches
- Titanium
- F9-003T (Titanium)



F9-011T Iris Hook and IOL Manipulator K-shape

- K-shape rotator tip
- 11mm angled shaft
- Overall length: 120mm, 4.7 inches
- Titanium



F9-013T Bechert Nucleus Rotator

- Y-shape rotator tip
- 11mm angled shaft
- Overall length: 120mm, 4.7 inches
- Titanium



F9-088T Kuglen Iris Hook

- Kuglen iris hook
- Overall length: 120mm, 4.7 inches
- Titanium



F9-121 Barrett Micro Mushroom Manipulator

- 11mm angled shaft 45 degree
- 0.4mm round tip
- Tip length 1.5mm
- Overall length: 120mm, 4.7 inches
- Stainless steel
- F9-121T (Titanium)



F9-121T Iris Retractors H-shape

- H-shape rotator tip
- 11mm angled shaft
- Overall length: 120mm, 4.7 inches
- Titanium



F9-123T Rentsch Boat Hook Straight

- Gamma-shape rotator tip
- 11mm angled shaft
- Overall length: 120mm, 4.7 inches
- Titanium



F9-151T Drysdale Nucleus Manipulator

- Flattened paddle tip
- Tip length 1.5mm
- Overall length: 120mm, 4.7 inches
- Titanium
- F9-151 (Stainless Steel)



F9-212T Anderson-Barrett Manipulator

- Tapered shafts
- Double ended instrument with chopper and barrett mushroom tip
- Overall length: 130mm, 5.1 inches
- Titanium



F9-218 Micro-Chop Lens Manipulator

- 90 degree curved and ball tip
- Overall length: 120mm, 4.7 inches
- Stainless steel



F9-222 Silverstein Phaco Manipulator & Quick Chop

- Universal slimmed profile crystalline lens
- Accommodate any phaco technique, or density of lens material
- Overall length: 120mm, 4.7 inches
- Stainless steel



F9-224 Sweeney Manipulator Cracker

- Angled for easy insertion into side port incision
- Flat, paddle for placement to crack the nucleus for emulsification
- Overall length: 114mm, 4.5 inches
- Stainless steel



F9-230T Lester Lens Manipulator Angled

- Angled shaft with 0.2mm hourglass shape tip
- Round knurled handle
- Overall length: 115mm, 4.5 inches
- Titanium



F9-232 Lindstrom Lens Positioner

- A semi-sharp tip to position an intraocular lens
- Overall length: 117mm, 4.6 inches
- Stainless Steel



F9-220 Osher Lens Manipulator

- Gently curved shaft
- Two prongs with desired chopping plane
- Overall length: 123mm, 4.9 inches
- Stainless steel



F9-301T Spatula Double End

- Double ended with 0.5mm by 12mm and 1mm by 10mm blades
- Both ends are blunt with gently curved blades
- Overall length: 128mm, 5 inches
- Titanium
- F9-301 (Stainless Steel)



F9-303T Iris Spatula 1.2mm x 10mm Angled

- Repositov
- 1.2mm x 10mm angled blade
- End 3.5 Blunt with gently curved blade
- Overall length: 118mm, 4.6 inches
- Titanium
- F9-303 (Stainless Steel)



F9-305 Iris Spatula 0.5mm x 10mm Angled

- Repositov
- 0.5mm x 10mm angled blade
- End 3.5 Blunt with gently curved blade
- Overall length: 118mm, 4.6 inches
- Stainless Steel
- F9-305T (Titanium)



F9-313 Spatula 1.2mm

- 1.2mm
- Overall length: 118mm, 4.6 inches
- Stainless Steel
- F9-315 (1.5mm)



F9-703 Kansas Nucleus Vectis Double End

- Double end for right and left handed surgeons
- 2.5mm x 12mm angled shafts
- Overall length: 130mm, 5.1 inches
- Stainless Steel
- F9-703T (Titanium)



F9-705 Kansas Nucleus Vectis Left / Right

- Left
- 2.5mmx12mm angled shaft
- 14mm from tip to bend
- Overall length: 120mm, 4.7 inches
- Stainless Steel
- F9-705-L (Left)
- F9-705-R (Right)



F9-711 Trabeculectomy Set

- 10mm long pointed tips with 3mm spread
- Overall length: 58mm, 2.3 inches
- Stainless Steel
- F9-711-L (Left)
- F9-711-R (Right)



F9-911T Flat Cover Double

- 8mm wide spoon
- Angled shaft
- Overall length: 113mm, 4.5 inches
- Titanium



F9-913T Lasik Flap Unzipper

- Double ended with 0.6mm by 5mm curved blades,
- 0.5mm by 6mm straight
- Overall length: 135mm, 5.3 inches
- Titanium



F9-929T Corneal Epithelium Spatula

- A paddle shaped tip which has broad
- Blunt edges on the opposite side which merge at the tip to create a penetrating point
- Overall length: 117mm, 4.6 inches
- Titanium



F9-935T Katzen Lasek Flap Unzipper

- Fan-shape tip
- 10mm angled shaft
- Overall length: 120mm, 4.7 inches
- Titanium



F9-604 Femto Lamellar / Smile Lamellar

- Dissector allows surgerys to open femto
- incisions prepared by the laser
- Overall length: 130mm, 5.1 inches
- Stainless Steel



F9-921T Corneal Epithelium Spatula

- A paddle shaped tip which has broad
- Blunt edges on one side and thin
- Sharp edges on the opposite side
- Overall length: 117mm, 4.6 inches
- Titanium
- F9-925T (Small)




F9-931T Femto Lasik Flap Spatula Double Ended

- 12mm curved shafts flat lifter and
- 1.25mm curved bullet shaped retreatment spatula
- Overall length: 127mm, 5 inches
- Titanium
- F9-931 (Stainless Steel)



F9-937T Katzen Lasek Flap Unzipper

- Angled shaft with 1.4mm wide and 4mm long quadrate tip
- 0.5mm tip angled
- Overall length: 120mm, 4.7 inches
- Titanium



F12-003 Smile Dissector

- Smile Double-Ended Dissector
- Overall length: 127mm, 5 inches
- Stainless Steel



F12-006

Smile Lenticule Splitter

- One end aims at helping surgeons to extract & remove the lenticule
- Overall length: 127mm, 5 inches
- Stainless Steel



Polisher

F9-605

Polisher Double-ended 27g

- 27g
- Designed for removing LECs from the inner surface of the capsule
- The circular ends are ideal for cleaning the equator of the capsular bag prior to lens insertion
- Provides a left angled and right angled polisher on one handle for surgeon convenience
- 1mm diameter both sides
- Overall length: 134mm, 5.3 inches
- Stainless Steel



I/A Handpiece and Flute Needle


Irrigating / Aspirating Handpiece

Flute Needle



F10-203T Bimanual I/A Handpiece Irrigation

- Dual port 20g
- Irrigating handpiece available in 0.8mm tip with front opening for minimal incision width
- Irrigating handpiece fits on male port of tubing
- Overall length: 120mm, 4.7 inches
- Titanium
- F10-203T-21 (21 gauge)
- F10-207T-23 (23 gauge)



F10-209T Bimanual I/A Handpiece Aspiration

- Single port 20g
- Aspiration handpiece fits on available in 0.8mm tip in curved models with a 0.3mm tip
- Overall length: 120mm, 4.7 inches
- Titanium
- F10-209T-21 (21 gauge)
- F10-211T-23 (23 gauge)



F10-213T I/A Handpiece Straight 20g Port

- 20g aspiration port
- Straight aspiration
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-213T-21 (21 gauge)
- F10-213T-23 (23 gauge)



F10-215T I/A Handpiece Angled 20g Port

- 20g aspiration port
- Angled aspiration
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-215T-21 (21 gauge)
- F10-215T-23 (23 gauge)



F10-217T I/A Handpiece Straight 45 degree 20g Port

- 20g aspiration port
- Straight 45 degree angled aspiration
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-217T-21 (21 gauge)
- F10-217T-23 (23 gauge)



F10-221T I/A Handpiece Straight for Mtp 20g Port

- 20g aspiration port
- I/A handpiece straight with silicone sleeve
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-221T-21 (21 gauge)
- F10-221T-23 (23 gauge)



F10-223T I/A Handpiece Angled for Mtp 20g Port

- 20g aspiration port
- I/A handpiece angled with silicone sleeve
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-223T-21 (21 gauge)
- F10-223T-23 (23 gauge)




F10-225T I/A Handpiece Straight 45 degree for Mtp 20g Port

- 20g aspiration port
- I/A handpiece straight 45 degree with silicone sleeve
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-225T-21 (21 gauge)
- F10-225T-23 (23 gauge)




F10-227T I/A Handpiece Straight for Alcon 20g Port

- 20g aspiration port
- Straight with silicone sleeve
- The blue irrigating sleeves on tip
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-227T-21 (21 gauge)
- F10-227T-23 (23 gauge)




F10-229T I/A Handpiece Angled for Alcon 20g Port

- 20g aspiration port
- Angled with silicone sleeve
- The blue irrigating sleeves on tip
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-229T-21 (21 gauge)
- F10-229T-23 (23 gauge)




F10-231T I/A Handpiece Straight 45 degree for Alcon 20g Port

- 20g aspiration port
- Straight 45 degree with silicone sleeve
- The blue irrigating sleeves on tip
- Overall length: 132mm, 5.2 inches
- Titanium
- F10-231T-21 (21 gauge)
- F10-231T-23 (23 gauge)





F10-401T Aspiration Flute Needle 20 gauge

- 20 gauge
- Uses 300mm silica cannulation
- Overall length: 158mm, 6.2 inches
- Titanium
- F10-403T (23 gauge)



F10-411 Aspiration Flute Needle

- Soft silicone brush tip mounted on 20 gauge
- Overall length: 35mm, 1.4 inches
- Stainless Steel
- F10-413 (23 gauge)



F10-421 Plastic tube

- Made of silica gel
- Overall length: 30mm, 1.2 inches
- Stainless Steel




11

Sterilizing Cases

Stainless Steel Tray

Sterilizing Tray



F11-001 Stainless Steel Tray

- Outside dimensions: 235mm x 148mm x 20mm
- Stainless steel tray for sterilizing instruments
- Silicone matting for supporting and protecting instruments
- F11-003 (215mm x 125mm x 18mm)




F11-005 Stainless Steel Tray Double Layer

- Outside dimensions: 235mm x 148mm x 35mm
- Stainless steel tray for sterilizing instruments
- Double silicone matting for supporting and protecting instruments




F11-021 Sterilizing Trays Double Layer

- Outside dimensions: 263mm x 162mm x 45mm
- Sterilizable plastic tray
- Double silicone matting for supporting and protecting instruments
- See-through lids




F11-023 Sterilizing Trays Large / Medium / Small

- Outside dimensions: 263mm x 162mm x 25mm
- Sterilizable plastic tray
- Silicone matting for supporting and protecting instruments
- See-through lids
- F11-025 (medium 198mm x 110mm x 20mm)
- F11-027 (small 158mm x 89mm x 25mm)




F11-041 Silicone Matting

- Dimensions: 247mm x 148mm x 16mm
- F11-043 (210mm x 130mm x 16mm)




F11-050 Sterilizing Trays

- Outside dimensions: 200mm x 110mm x 45mm
- Sterilizable plastic tray
- Silicone matting for supporting and protecting instruments
- See-through lids




F11-052 Sterilizing Trays

- Outside dimensions: 200mm x 110mm x 45mm
- Sterilizable plastic tray
- For 3 phaco handpieces with special silicone plug-in inserts
- See-through lids



F11-054 Sterilizing Trays

- Outside dimensions: 265mm x 162mm x 45mm
- Sterilizable plastic tray
- For 3 phaco handpieces with special silicone plug-in inserts
- See-through lids



Trypan Blue

AlgerBrush Burr
Trypan Blue
Iris Retractors



F12-AB Trypan Blue Auroblue

- Trypan blue solution 0.06% W/v - 1ml
- 5 vial / box



F12-BB Trypan Blue Bioblue

- Trypan blue ophthalmic solution 0.06% W/v
- Use as an adjunct in ophthalmic surgery
- 10 vial / box
- 1ml / vial



F12-BC Trypan Blue Biochol

- Carbachol
- 10 vial / box
- 1ml / vial




F12-043 AlgerBrush Burr

- Black flat occluder 63mm x 52mm
- Round knurled handle
- Overall length 200mm, 7.9 inches
- 0.5mm / 1.0mm / 2.5mm / 3.5mm / 4.0mm



F12-119 Reusable Iris Retractors


- 10mm long
- 1 x 6 pcs
- Bright blue micro hooks




Blade Breaker

Disposable Bard Parker

Sapphire Knife



Sapphire Knife



F13-021 Sapphire Knife Side Port

- 15 degree
- 1.0mm
- Titanium Handle




F13-023 Sapphire Knife Spear

- C-45
- 1.0mm
- Titanium Handle




F13-025 Sapphire Knife Lancet 3.2mm

- E-1
- 3.2mm
- Titanium Handle
- F13-027 (3.0mm)
- F13-033 (2.75mm)
- F13-035 (2.2mm)
- F11-037 (2.8mm)




F13-029 Sapphire Knife Crescent

- E-2
- Titanium Handle



F13-031 Sapphire Knife Trapezoid

- E-3
- Titanium Handle



F13-505T Swiss Model Blade Breaker Straight 10mm Serrated Jaws

- Features conical 10mm serrated jaws
- Round knurled handle with enf lock
- Overall length: 117mm, 4.6 inches
- Titanium
- F13-505 (stainless steel)



F13-507 Swiss Model Blade Breaker Angled 9mm Serrated Jaws

- Blade breaker and holder
- Angle 9mm serrated jaws
- Overall length: 118mm, 4.6 inches
- Stainless Steel



F13-511 Bard Parker Handle

- For disposable slot blades
- 3#
- Stainless Steel



F13-519 Disposable Bard Parker Round Tip

- Round tip
- Sterile
- 15#
- Stainless Steel
- F13-515 (sharp tip)



Surgical Package Set

Cataract Surgery Set

Phacoemulsification Set

Glaucoma Set

Strabismus Surgery Set

Lasek Surgery Set

Blepharoplasty Set

Dacryocystorhinostomy Set




Stainless Steel

1. F1-003	Wire Speculum	p.7	12. F3-003	Bonn Suturing Forceps Straight	p.15
2. F2-081	Infusion Cannula Angled 18g	p.12	13. F3-105	McPherson Tying Forceps Straight	p.22
3. F2-113	Irrigation/ Aspiration Cannula 22g	p.12	14. F3-109	McPherson Tying Forceps Angled	p.22
4. FO62119	Coaxial I/A Cannula	p.---	15. F3-303	Intraocular Lens Forceps Curved	p.28
5. F2-091	Infusion Lens Loop 18g	p.12	16. F3-401	Pierse Notched Forceps	p.19
6. F9-001	Sinskeys Hook	p.59	17. F4-031	Vannas Scissors Straight	p.36
7. F9-003	Fenzl Hook	p.59	18. F4-033	Vannas Scissors Curved	p.36
8. FO69831	Cautery	p.---	19. F4-001	Corneal Scissors Curved	p.34
9. F9-301	Specula Double Ended	p.61	20. F4-047	Vannas Csapsulotomy Scissors	p.36
10. FO62085	Infusion Handpiece	p.---	21. F5-001	Barraquer Needle Holder	p.40
11. F9-505	Swiss Model Blade Breaker Straight	p.74			

Titanium

1. F1-003T	Wire Speculum	p.7	12. F3-003T	Bonn Suturing Forceps Straight	p.15
2. F2-081T	Infusion Cannula Angled 18g	p.12	13. F3-105T	McPherson Tying Forceps Straight	p.22
3. F2-113T	Irrigation/ Aspiration Cannula 22g	p.12	14. F3-109T	McPherson Tying Forceps Angled	p.22
4. FO62119	Coaxial I/A Cannula	p.---	15. F3-303T	Intraocular Lens Forceps Curved	p.28
5. F2-091T	Infusion Lens Loop 18g	p.12	16. F3-401T	Pierse Notched Forceps	p.19
6. F9-001T	Sinskeys Hook	p.59	17. F4-031T	Vannas Scissors Straight	p.36
7. F9-003T	Fenzl Hook	p.59	18. F4-033T	Vannas Scissors Curved	p.36
8. FO69831	Cautery	p.---	19. F4-001T	Corneal Scissors Curved	p.34
9. F9-301T	Specula Double Ended	p.61	20. F4-047T	Vannas Csapsulotomy Scissors	p.36
10. FO62085	Infusion Handpiece	p.---	21. F5-001T	Barraquer Needle Holder	p.40
11. F9-505T	Swiss Model Blade Breaker Straight	p.74			




Stainless Steel

1. F1-005	Wire Speculum	p.7	8.F9-103	Nagahra Nucleus Chopper 45°	p.56
2. F11-021	Sapphire Knife Side Port	p.73	9.F9-113	Nagahra Nucleus Chopper 90°	p.57
3. F11-029	Sapphire Knife Crescent	p.73	10.F3-217	Capsulorhexis Forceps	p.20
4. F11-025	Sapphire Knife Lancet	p.73	11.F3-323	Intraocular Lens Folding Forceps	p.29
5. F9-011	Iris Hook and IOL Manipulator	p.60	12.F3-327	Lens Implantation Forceps	p.29
6. F9-013	Becherf Nucleus Rotator	p.60	13.F10-213	I/A Handpiece Straight 20g Port	p.65
7. F9-088	Kuglen Iris Hook	p.60	14.F10-221	I/A Handpiece Straight for Mtp 20g Port	p.65

Titanium

1. F1-005T	Wire Speculum	p.7	8.F9-103T	Nagahra Nucleus Chopper 45°	p.56
2. F11-021T	Sapphire Knife Side Port	p.74	9.F9-113T	Nagahra Nucleus Chopper 90°	p.57
3. F11-029T	Sapphire Knife Crescent	p.74	10.F3-217T	Capsulorhexis Forceps	p.20
4. F11-025T	Sapphire Knife Lancet	p.74	11.F3-323T	Intraocular Lens Folding Forceps	p.29
5. F9-011T	Iris Hook and IOL Manipulator	p.60	12.F3-327T	Lens Implantation Forceps	p.29
6. F9-013T	Becherf Nucleus Rotator	p.60	13.F10-213T	I/A Handpiece Straight 20g Port	p.65
7. F9-088T	Kuglen Iris Hook	p.60	14.F10-221T	I/A Handpiece Straight for Mtp 20g Port	p.65



Stainless Steel


1. F6-011	Castroviejo Caliper Straight Shafts	p.43
2. F4-076	Eye Scissors Curved	p.37
3. F4-111	Baeaaquer Iris Scissors	p.37
4. F4-101	Corneal Scissors	p.35
5. F3-003	Bonn Suturing Forceps Straight	p.15
6. F3-551	Panfique Suturing Forceps Straight	p.19
7. F7-023	Descemet Membrane Punch	p.49
8. F9-711-L	Trabeculectomy Set (left)	p.61
9. F9-711-R	Trabeculectomy Set (right)	p.61
10. F9-301	Spatula Double Ended	p.61
11. F9-921	Corneal Epithelium Spatula	p.62
12. F9-315	Spatula 1.5mm	p.61
13. FO69831	Cautery	p.----

Titanium

1. F6-011T	Castroviejo Caliper Straight Shafts	p.43
2. F4-076T	Eye Scissors Curved	p.37
3. F4-111T	Baeaaquer Iris Scissors	p.37
4. F4-101T	Corneal Scissors	p.35
5. F3-003T	Bonn Suturing Forceps Straight	p.15
6. F3-551T	Panfique Suturing Forceps Straight	p.19
7. F7-023T	Descemet Membrane Punch	p.49
8. F9-711T-L	Trabeculectomy Set (left)	p.61
9. F9-711T-R	Trabeculectomy Set (right)	p.61
10. F9-301T	Spatula Double Ended	p.61
11. F9-921T	Corneal Epithelium Spatula	p.62
12. F9-315T	Spatula 1.5mm	p.61
13. FO69831	Cautery	p.----




1. FO81013	Kershner Eye Speculum	p.---
2. F1-015T	Kershner Eye Speculum	p.8
3. F3-105	McPherson Tying Forceps Straight 5mm	p.22
4. F3-256	Castroviejo Suturing Forceps Straight	p.16
5. F3-505-L	Muscle Clamp Curved (left)	p.31
6. F3-505-R	Muscle Clamp Curved (right)	p.31
7. F3-509	Muscle Clamp Straight	p.31
8. F4-101	Corneal Scissors	p.36
9. F4-053	Tenotomy Scissors Curved 20mm	p.37
10. F5-017	Castroviejo Needle Holder with Lock	p.41
11. F8-001	Desmarres Lid Retractor 8mm	p.53
12. F8-093	Jameson Strabismus Hook 1.5mm	p.52
13. F8-091	Jameson Strabismus Hook 2.5mm	p.52
14. F8-089	Jameson Strabismus Hook 6mm	p.52
15. F8-085	Jameson Strabismus Hook 12mm	p.52
16. F8-083	Jameson Strabismus Hook 10mm	p.52
17. F8-103	Lid Retractor 7.2mm	p.53
18. F8-105	Lid Retractor 9.2mm	p.53
19. FO68099	Jameson Strabismus Hook Double arm	p.---
20. F8-095	Jameson Strabismus Hook 0.7mm	p.53
21. F8-097	Jameson Strabismus Hook 1.75mm	p.52
22. F8-061	Jameson Strabismus Hook 6mm	p.52
23. F8-069	Jameson Strabismus Hook 6mm with hole	p.52
24. F8-077	Jameson Strabismus Hook 8mm	p.52



Titanium

1. FO62013	Vidaurri Lasek Cannula	p.---
2. F9-902T	Outside Optic Zone Marker 8mm	p.45
3. F9-905T	Inside optic Zone Marker 8.5mm	p.46
4. F9-921T	Corneal Epithelium Spatula	p.62
5. F9-925T	Small Corneal Epithelium Spatula	p.62
6. F9-929T	Corneal Epithelium Spatula	p.62
7. F9-935T	Katzen Lasek Flap Unzipper	p.62
8. F9-937T	Katzen Lasek Flap Unzipper	p.62
9. F9-931T	Femto Lasik Flap Spatula Double Ended	p.62
10. F9-915T	Castroviejo Spatula	p.---
11. F9-917T	Castroviejo Spatula Double End	p.---
12. F9-913T	Lasik Flap Unzipper	p.62
13. F9-911T	Flat Cover Spoon	p.61



1. F9-505	Swiss Model Blade Breaker Straight	p.74
2. F6-011T	Castroviejo Caliper Straight Shafts	p.41
3. F5-063	Jones Towel Clamp	p.41
4. F3-005-90	Castroviejo Suturing Forceps 0.12mm Straight	p.16
5. F3-551T	Paufique Suturing Forceps Straight	p.19
6. F3-553T	Paufique Suturing Forceps Curved	p.19
7. F3-715	Suturing Bishop-Harmon Tissue Forceps	p.19
8. F5-039	Halstead Straight Hemostatic Mosquito Forceps	p.41
9. F5-057	Halstead Curved Hemostatic Mosquito Forceps	p.41
10. F3-588T	Putterman clamp	p.30
11. F8-061	Jameson Hook 6mm	p.52
12. F9-849	Jaeger Lid Plate	p.45
13. F5-007T	Barraquer Needle Holder with Lock	p.40
14. F3-761	Fixation Graefe Fixation Forceps	p.27
15. F8-001	Desmarres Lid Retractor 8mm	p.53
16. F8-003	Desmarres Lid Retractor 10mm	p.53
17. F8-005	Desmarres Lid Retractor 12mm	p.53
18. F4-053	Tenotomy Scissors Curved 20mm	p.37
19. F4-057	Bonn Model Eye Scissors	p.37
20. F4-085	Steven Straight Tenotomy Scissors	p.38
21. F10-023	Sterilizing Trays Large	p.68



1. F7-007	Kerrison Rongeur 2mm x 3mm 100 degree	p.48
2. FO61081	Lacrimal Sac Retractor	p.----
3. F9-505	Blade Breaker	p.74
4. F3-701	Dressing Forceps Straight	p.26
5. F3-256-45	Castroviejo Suturing Forceps 0.3mm Straight	p.16
6. F5-021	Castroviejo Needle Holder 10mm	p.41
7. F4-072	Eye Scissors Straight	p.37
8. F5-039	Halstead Straight Hemostatic Mosquito Forceps	p.41
9. F5-057	Halstead Curved Hemostatic Mosquito Forceps	p.27
10. F2-143	Lacrimal Probes	p.13
11. F2-171	Lacrimal Dilators	p.54
12. FO68031	Fixation Hook	p.----
13. FO69865	Freer Lacrimal Chisel (Big)	p.----
14. FO69867	Freer Lacrimal Chisel (Small)	p.----




About Us

Our factory has been established for over 20 years, with the vision of providing the surgical world with highly skilled and advanced manufacturing of surgical instruments.

We consult with surgeons from all around the world and develop innovative instruments for new surgical techniques and advances in the field of ophthalmology. Combined our handcraft tradition with technological developments and equipped with high precision CNC and automation machine, our qualified working teams guarantee the consistency of the products quality.

All our instruments are CE registered and ISO13485:2016 certified.



Your Care, Our Focus



For the most up-to-date listing of available instruments, visit:
www.focusinstruments.com.hk

Network Technology
Development Ltd.

CE MARK
FDA APPROVED
ISO13485 CERTIFIED

Unit 21, 10/F, Yale Ind Center,
61-63 Au Pui Wan Street,
Fotan, Hong Kong



Shop Online

Customer Service:

+852-3590 5727
Fax: +852-3012 1635
Email: alfredau0310@hotmail.com
zhengchfr@hotmail.com