


Facility Medical Co.


Suture List




SURGICRYL
PGA
POLYGLYCOLIC ACID




POLYPROPYLENE
BLUE



DACLON
NYLON
BLACK



SILK
BRAIDED
BLACK



SILK
VIRGIN
BLUE



POLYESTER



MEDICAL DEVICES & OPHTHALMOLOGY TRADING

📞 +966 54 554 9351 - +966 53 438 3566

📍 12891, Riyadh, Kingdom of Saudi Arabia

🌐 <https://facilitymedical.com>

✉ sales@facilitymedical.com

ITEM LIST



Facility Medical Co.
Nadah Rd, Riyadh, Saudi Arabia
+966 54 554 9351 +966 53 438 3566
sales@facilitymedical.com
<https://facilitymedical.com>

Surgical Suture

Product Code	Product Name	Description	Product Image	Origin	Price
653	PGA <i>Rev. Cutting</i>	SURGICRYL PGA EP 1 USP 5/0 NEEDLE 3/8 CIR REV. CUTTING 12MM-45cm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
672	PGA <i>Circ. Spatulated</i>	SURGICRYL PGA EP1 USP 5/0 45cm 2 NEEDLES 1/4 CIR SPATULATED 8MM 350µm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
651	PGA <i>Cir Taper</i>	SURGICRYL PGA EP 0.7 USP 6/0 NEEDLE 1/2 CIR TAPER 13MM-75cm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
652	PGA <i>Rev. Cutting</i>	SURGICRYL PGA EP 0.7 USP 6/0 NEEDLE 3/8 CIR REV. CUTTING 12MM- 45cm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
674	PGA <i>Circ. Spatulated</i>	SURGICRYL PGA EP0.7 USP 6/0 45cm 2NEEDLES 1/4CIR SPATULA 8MM 300µm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
675	PGA <i>Circ. Spatulated</i>	SURGICRYL PGA EP0.4 USP 8/0 45cm 2NEEDLES 3/8CIR SPATULA 6MM 200µm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
676	PGA <i>Circ. Spatulated</i>	SURGICRYL PGA EP 0.3 USP 9/0 30cm 2 NEEDLES 3/8CIR SPATULA 6.4 MM 200µm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
677	PGA <i>Circ. Spatulated</i>	SURGICRYL PGA EP 0.2 USP 10/0 30cm 2 NEEDLES 3/8CIR SPATULA 6.2MM 150µm - <u>ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
657	POLYPROPYLENE <i>Cir Taper</i>	POLYPROPYLENE BLUE 75cm EP 1 USP 5/0 DOUBLE NEEDLE 1/2 CIR TAPER 13 <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
669	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 6/0 45cm 2 X 1/4 CIRC.SPATULED 8 MM 300µm - <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request

ITEM LIST

Surgical Suture

Product Code	Product Name	Description	Product Image	Origin	Price
656	POLYPROPYLENE <i>Cir Taper</i>	POLYPROPYLENE BLUE 60cm EP 0.7 USP 6/0 2 NEEDLES 1/2 CIR TAPER 13 MM - <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
655	POLYPROPYLENE <i>Cir Taper</i>	POLYPROPYLENE BLUE 60cm EP0.5 USP 7/0 2 x 3/8 CIR TAPER POINT 9.3 MM <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
679	POLYPROPYLENE <i>Cir Taper</i>	POLYPROPYLENE BLUE 60cm EP 0.5 USP 7/0 2 NEEDLES 1/2 CIR TAPER 8MM <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
654	POLYPROPYLENE <i>Cir Soft Coroner</i>	POLYPROPYLENE BLUE 60cm EP0.4 USP 8/0 2 x 3/8 CIR SOFT CORONAR PT 8MM <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
668	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 9/0 30cm 2 NEEDLES 3/8 CIRC.SPATULED 6.40 MM 200µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
668	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 9/0 30cm 2 NEEDLES 3/8 CIRC.SPATULED 6.40 MM 200µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
665	POLYPROPYLENE <i>Straight Spatulated</i>	POLYPROPYLENE USP 10/0 20cm 2STRAIGHT NEEDLES 16.15 MM 150µm SPATULATED <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
666	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 10/0 20cm 2 X 1/4 CIRC.SPATULED 16 MM 150µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
666	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 10/0 20cm 2 X 1/4 CIRC.SPATULED 16 MM 150µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
667	POLYPROPYLENE <i>Circ. Spatulated</i>	POLYPROPYLENE USP 10/0 30cm 2 X 3/8 CIRC.SPATULED 6.40 MM 200µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request

ITEM LIST

Surgical Suture

Product Code	Product Name	Description	Product Image	Origin	Price
658	NYLON <i>CIR REV CUT.</i>	DACLON NYLON BLACK 75cm EP 0.7 USP 6/0 NEEDLE 3/8 CIR REV CUT. 12 MM <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
662	NYLON <i>Circ. Spatulated</i>	DACLON NYLON BLACK USP 7/0 30cm 2 X 3/8 CIRC.SPATULED 6.20 MM 200µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
663	NYLON <i>Circ. Spatulated</i>	DACLON NYLON BLACK USP 8/0 30cm 2 NEEDLES 3/8 CIRC.SPATULED 6.20 MM 150µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
664	NYLON <i>Circ. Spatulated</i>	DACLON NYLON BLACK USP 9/0 30cm 2 NEEDLES 3/8 CIRC.SPATULED 6.20 MM 150µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
660	NYLON <i>Circ. Spatulated</i>	DACLON NYLON BLACK USP 10/0 30cm 2 X 3/8 CIRC.SPATULED 6.40 MM 220µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
661	NYLON <i>Circ. Spatulated</i>	DACLON NYLON BLACK USP 10/0 38cm 2 X 3/8 CIRC.SPATULED 6.40 MM 220µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
670	SILK <i>Circ. Spatulated</i>	SILK BLACK USP 6/0 45cm 2 X 3/8 CIRC.SPATULED 6.5 MM 240µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
671	SILK <i>Circ. Spatulated</i>	SILK VIRGIN BLUE USP 8/0 45cm 2 NEEDLES 3/8 CIRC.SPATULED 6.5 MM 200µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request
672	POLYSTER <i>Circ. Spatulated</i>	POLYESTER USP 5/0 45cm 2 NEEDLES 1/4 CIRC.SPATULED 8 MM 350µm <u>NON-ABSORBABLE</u> - 12 Sutures		Belgium	As Per Request

ITEM LIST

S u r g i c a l S u t u r e

Product Code	Product Name	Description	Product Image	Origin	Price
--------------	--------------	-------------	---------------	--------	-------

SMI surgical needles are manufactured from highest quality stainless steel grades for outstanding performance, combining high ductility and strong bending moment. Advanced silicon coating ensures smooth penetration and reduced tissue trauma. The shape of the needle body is designed in order to allow maximum stability in the needle holder.

5/8 CIRCLE



1/2 CIRCLE



3/8 CIRCLE




1/4 CIRCLE



J-NEEDLE




STRAIGHT




TAPER POINT

A gradual taper to a round bodied needle.




SIDE CUTTING SPATULA NEEDLES

Needle with cutting edges on the two sides.




TAPER CUTTING POINT

Balanced, cutting point-tapering into a round bodied needle.




SOFT CUTTING POINT - CORONARY

A minute point followed by a round body and a rectangular main body for better penetration in cardiovascular surgery.




REVERSE CUTTING

A triangular body offers a third cutting edge on the outer convex curvature.




CONVENTIONAL CUTTING

A triangular body offers a third cutting edge on the inside concave curvature.



EASY PASS NEEDLE

Precision point with extra sharpening.



BLUNT NEEDLE

Cylindrical bodied needle with a rounded point. For soft tissue such as in liver and kidneys. Reduced risk of needle stick injuries.